

 **ENHANCE**[™]
TRANSCAROTID ACCESS KIT
PERIPHERAL


Marks of Confidence Inside and Out



SILKROAD
MEDICAL[®]

SUPPORTIVE 0.018" MARKED WIRE

- Distal support for smoother microsheath entry
- Pre-shaped tip for selective navigation
- Marked wire for precise placement in vasculature



Wire mark
Indicates 3-5 cm
of wire in vessel

5cm

3cm



4F MARKED MICROSHEATH

- Smooth dilator to sheath transition
- Depth markings for accurate positioning
- Pre-loaded with non-stiffened dilator



EXTENSION TUBE

- 20cm extension tube for contrast injections
- Complete with 3-way stopcock



Needle mark
5mm from beveled tip


Pre-shaped tip



“The TCAR procedure hinges on arterial access.”

Richard Cambria, MD

A Complete Solution for TCAR



ENHANCE™ Transcarotid Peripheral Access Kit

Catalog Number	Components	Shipped
SR-4F21G7D-MP	7cm 21G needle 50cm 0.018" microwire 4F 15cm sheath & dilator 4F stiffened dilator 20cm extension tube	10/box

Manufactured by Galt Medical for Silk Road Medical distribution.

Indications for Use: These access kits are intended to introduce up to a 0.038" guidewire or catheter into the peripheral vascular system, including the carotid artery, following a small gauge needle stick.

Caution: Federal (U.S.) law restricts this device to sale by or on the order of a physician.



SILK ROAD MEDICAL, INC.
1213 INNSBRUCK DRIVE
SUNNYVALE, CA 94089

TOLL-FREE: 855.410.TCAR (8227)
MAIN: 408.720.9002
FAX: 408.720.9013